



CareOne MediMart



## G3 BiPAP System B25VT

A bi-level device with S/T mode and target tidal volume function for patients with respiratory insufficiency.

Working mode & functions: CPAP, S, T, S/T, Target Tidal Volume  
Pressure range: 4 - 25 cmH<sub>2</sub>O  
Advanced event detection: CSA Detection  
Auto Ramp: Yes  
Auto Humidifier: Yes  
Auto Heated Tubing: Yes  
Reslex (Expiration Pressure Release): Yes

- Auto Leak Compensation
- Inspiratory time control
- PM2.5 filter box
- A 3.5" screen displays real-time data
- Cellular / WiFi Module

Various ways of reviewing Sleep Report like Device screen, PAP Link App, Cloud Platform, iCode web version ([www.bmc-icode.com](http://www.bmc-icode.com)), RESmart nPAP Data Analysis Software (PC software) RESmart Software web version ([www.icodeconnect.com/quick/info](http://www.icodeconnect.com/quick/info))

BMC Authorised Service Centre





 CONTACT US

### BMC Medical Co. Ltd.

Room 110, Tower A Fengyu Building, No.115 Fucheng Road,  
Haidian, 100036 Beijing, PEOPLE'S REPUBLIC OF CHINA  
Tel: +86-10-5166 3880  
Fax: +86-10-5166 3880 ext. 810  
[en.bmc-medical.com](http://en.bmc-medical.com)  
SC-02-21-06-004-Ver2.1

### Bhupinder Nanda

A-4/1, Sector 18, Rohini, Delhi - 110089, INDIA

 +918860858842  +918459072606

 [shopcareone.com](http://shopcareone.com)

 [careonedelhi@gmail.com](mailto:careonedelhi@gmail.com)